



**Institute of
Applied Physics**

Friedrich-Schiller-Universität Jena

Microscopy

Lecture 2: Optical System of the Microscopy II

2012-10-22

Herbert Gross

2 Optical System of the Microscopy II

Preliminary time schedule

No	Date	Main subject	Detailed topics	Lecturer
1	15.10.	Optical system of a microscope I	overview, general setup, binoculars, objective lenses, performance and types of lenses, tube optics	Gross
2	22.10.	Optical system of a microscope II	Etendue, pupil, telecentricity, confocal systems, illumination setups, Köhler principle, fluorescence systems and TIRF, adjustment of objective lenses	Gross
3	29.10.	Physical optics of widefield microscopes	Point spread function, high-NA-effects, apodization, defocussing, index mismatch, coherence, partial coherent imaging	Gross
4	05.11.	Performance assessment	Wave aberrations and Zernikes, Strehl ratio, point resolution, sine condition, optical transfer function, conoscopic observation, isoplanatism, straylight and ghost images, thermal degradation, measuring of system quality	Gross
5	12.11.	Fourier optical description	basic concepts, 2-point-resolution (Rayleigh, Sparrow), Frequency-based resolution (Abbe), CTF and Born Approximation	Heintzmann
6	19.11.	Methods, DIC	Rytov approximation, a comment on holography, Ptychography, DIC	Heintzmann
7	26.11.	Imaging of scatter	Multibeam illumination, Cofocal coherent, Incoherent processes (Fluorescence, Raman), OTF for incoherent light, Missing cone problem, imaging of a fluorescent plane, incoherent confocal OTF/PSF	Heintzmann
8	03.12.	Incoherent emission to improve resolution	Fluorescence, Structured illumination, Image based identification of experimental parameters, image reconstruction	Heintzmann
9	10.12.	The quantum world in microscopy	Photons, Poisson distribution, squeezed light, antibunching, Ghost imaging	Wicker
10	17.12.	Deconvolution	Building a forward model and inverting it based on statistics	Wicker
11	07.01.	Nonlinear sample response	STED, NLSIM, Rabi the information view	Wicker
12	14.01.	Nonlinear microscopy	two-photon cross sections, pulsed excitation, propagation of ultrashort pulses, (image formation in 3D), nonlinear scattering, SHG/THG - symmetry properties	Heisterkamp
13	21.01.	Raman-CARS microscopy	principle, origin of CARS signale, four wave mixing, phase matching conditions, epi/forward CARS, SRS.	Heisterkamp
14	28.01	Tissue optics and imaging	Tissue optics, scattering&aberrations, optical clearing,Optical tomography, light-sheet/ultramicroscopy	Heisterkamp
15	04.02.	Optical coherence tomography	principle, interferometry, time-domain, frequency domain.	Heisterkamp

2 Optical System of the Microscopy II

Contents of 2nd Lecture





1. Etendue
2. Pupil
3. Telecentricity
4. Confocal systems
5. Illumination setups
6. Köhler principle
7. Fluorescence systems and TIRF
8. Objective adjustment

2 Optical System of the Microscopy II

Microscope Objective Lens

- No rigid relationship between magnification and aperture
- Product of field size and NA fixes the overall capacity
- Variation of image-sided numerical aperture




2 Optical System of the Microscopy II

Etendue of Microscope Objective Lenses

- Information capacity, etendue, space-bandwidth product, number of PSF's resolved in the field

$$G = \frac{\pi}{4} \cdot (D_{field} \cdot NA)^2$$

- Non-rigid correlation between magnification and etendue
- Largest etendue for medium magnifications
- Immersion systems have larger etendue



2 Optical System of the Microscopy II


Properties of the Pupil

Relevance of the system pupil :

- Brightness of the image
Transfer of energy
- Resolution of details,
Information transfer, location of Fourier spectrum
- Image quality
Aberrations due to aperture
- Image perspective
Perception of depth
- Compound systems:
matching of pupils is necessary, location and size

2 Optical System of the Microscopy II



Entrance and Exit Pupil



2 Optical System of the Microscopy II

Microscope Objective Lens: Pupil

- Object space telecentric
- Real rear stop is not defining the pupil
- Collimated outgoing beam
- Exit pupil usually not accessible




2 Optical System of the Microscopy II





Microscope Objective Lens: Pupil

- Diameter of pupil only weak correlated with magnification
- Pupil distortion:

$$D_{distExp} = \frac{y'_p}{f' \cdot \sin u} - 1$$





- Pupil distortion small (sine condition corrected)



2 Optical System of the Microscopy II

Microscope Objective Lens: Pupil



- Imaging of the pupil is important
- Residual aberrations : sharp edge of aperture is desired
- Real systems : smooth illumination profile



2 Optical System of the Microscopy II

Microscope Objective Lens: Pupil


- Example of larger residual aberrations of pupil image :
 - axial shift of pupil with field size
 - no sharp imaging of aperture stop
 - coloured edge of stop boundary



2 Optical System of the Microscopy II

Microscope Objective Lens: Pupil



- Vignetting of the pupil :
 1. truncation of the bundle for finite field sizes
 2. chief ray not identical with centroid
 3. perturbation of telecentricity
 4. in microscopic systems mostly at the front lens



2 Optical System of the Microscopy II

Vignetting



- Truncation at front lens
- Pupil vignetting:
crescent shape of
light cone
- Illumination fall-off towards
field boundary



2 Optical System of the Microscopy II

Vignetting


- Truncation of regions with large aberrations as correction method
- Improved performance
- Psf elliptical, anisotropic resolution
- Energy reduced



2 Optical System of the Microscopy II

Telecentricity




- Special stop positions:
 1. stop in back focal plane: object sided telecentricity
 2. stop in front focal plane: image sided telecentricity
 3. stop in intermediate focal plane: both-sided telecentricity
- Telecentricity:
 1. pupil in infinity
 2. chief ray parallel to the optical axis



2 Optical System of the Microscopy II

Telecentricity


- Design problems with telecentricity:
Usual telecentricity only fulfilled for one wavelength
- Telecentricity is related to the centroid ray.
Therefore the telecentricity is disturbed by vignetting effects



2 Optical System of the Microscopy II

Confocal Microscope



- Laser scan microscope
- Depth resolution (sectioning) with confocal pinhole
- Transverse scan on field of view
Digital image
- Only light coming out of the conjugate plane is detected
- Perfect system: scan mirrors conjugate to pupil location
- System needs a good correction of the objective lens, symmetric 3D distribution of intensity



2 Optical System of the Microscopy II

Confocal Images


Depth resolved images



2 Optical System of the Microscopy II




Confocal Laser Scan Microscope

- Complete setup: objective / tube lens / scan lens / pinhole lens
- Scanning of illumination / descanning of signal
- Scan mirror conjugate to system pupil plane
- Digital image processing necessary



2 Optical System of the Microscopy II Confocal Laser Scan Microscope



- Scan lens
- Diffraction limited
- Change in pupil location of objective lenses is critical perturbation of telecentricity



2 Optical System of the Microscopy II

Confocal Laser Scan Microscope

- Pinhole lens
- Only axial colour is essential
- Usage on axis only due to descanning
- Variable pinhole size not too small:
small aperture, retrofocus lens



2 Optical System of the Microscopy II

Confocal Microscopy: PSF and Lateral Resolution

- Normalized transverse coordinate v

$$v = \frac{2\pi}{\lambda} \cdot x' \cdot \sin \alpha$$


- Usual PSF: Airy

$$I(v) = \left[\frac{2J_1(v)}{v} \right]^2$$

- Confocal imaging:
Identical PSF for illumination and observation
assumed

$$I(v) = \left[\frac{2J_1(v)}{v} \right]^4$$

Resolution improvement by factor 1.4 for FWHM



2 Optical System of the Microscopy II

Confocal Microscopy: Axial Sectioning

- Normalized axial coordinate

$$u = \frac{8\pi}{\lambda} \cdot z \cdot \sin^2(\alpha/2)$$

- Conventional wide field imaging:

Intensity on axis

$$I(u) = \left[\frac{\sin(u/2)}{u/2} \right]^2$$

Axial resolution

$$\Delta z_{wide}^{(approx)} = \frac{0.45 \cdot \lambda}{n' \cdot (1 - \cos \theta)}$$


- Confocal imaging:

Intensity on axis

$$I(u) = \left[\frac{\sin(u/2)}{u/2} \right]^4$$

Axial resolution improved by factor 1.41 for FWHM

$$\Delta z_{confo} = \frac{0.319 \cdot \lambda}{n' \cdot (1 - \cos \theta)}$$



2 Optical System of the Microscopy II

Size of Pinhole and Confocality



- Large pinhole: geometrical optic
- Small pinhole:
 - Diffraction dominates
 - Scaling by Airy diameter $a = D/D_{\text{Airy}}$
 - diffraction relevant for pinholes

$$D < D_{\text{airy}}$$

- Confocal signal:

Integral over pinhole size


$$S(u) = \int_0^a |U(u, v)|^2 2\pi v dv$$



2 Optical System of the Microscopy II


Confocal Signal with Spherical Aberration

- Spherical aberration:
 - PSF broadened
 - PSF no longer symmetrical around image plane during defocus
- Confocal signal:
 - loss in contrast
 - decreased resolution



2 Optical System of the Microscopy II Confocal Laser Scan Microscope


- Depth signal as a function of wavelength
- Disturbance by axial colour aberration



2 Optical System of the Microscopy II

Illumination Optics: Overview

- Four possibilities for practical needs
- Epi vs. trans-illumination
- Bright vs. dark field illumination
- Comparison of light cones for imaging and illumination parts




2 Optical System of the Microscopy II


Illumination Optics: Overview

- Instrumental realizations


a) incident illumination bright field




b) incident illumination dark field



c) transmitted illumination bright field



d) transmitted illumination dark field




2 Optical System of the Microscopy II

Illumination Optics: Overview

- Typical images for different illuminations

epi



bright



trans



dark





2 Optical System of the Microscopy II

Köhler Illumination Principle

Principle of Köhler illumination:

- Alternating beam paths of field and pupil
- No source structure in image
- Light source conjugated to system pupil
- Differences between ideal and real ray paths



2 Optical System of the Microscopy II

Illumination Optics: Overview

- Types of settings :

1. Köhler :

source into pupil,
mostly used


2. Critical :

source into field of view,
source structure disturbs image

3. Projection :

source into condenser




4. Arbitrary



2 Optical System of the Microscopy II

Illumination Optics: Collector


- Requirements and aspects:
 1. Large collecting solid angle
 2. Correction not critical
 3. Thermal loading large
 4. Mostly shell-structure for high NA




2 Optical System of the Microscopy II

Illumination Optics: Condenser

2. Abbe type, achromatic, NA = 0.9 , aplanatic, residual spherical



3. Aplanatic achromatic, NA = 0.85




2 Optical System of the Microscopy II

Illumination Optics: Condenser



- Dark-field illumination systems

- 1. Trans-illumination

- Cardioid-mirror



Realizations: approximation of Cardioid curve





2 Optical System of the Microscopy II

Illumination Optics: Condenser

2. Epi-illumination



Complicated ring-shaped components around objective lens



2 Optical System of the Microscopy II

Illumination Optics: TIRF-Illumination


- Epi- and trans illumination for TIRF




2 Optical System of the Microscopy II

Fluorescence Microscopy

- Fluorescence microscopy is the most frequently employed mode of light microscopy used in biomedical research today
- Setup:




- Necessary components:
Dichroic beam splitter, excitation filter with sharp edge




2 Optical System of the Microscopy II

TIRF Microscopy



- Total internal reflection microscopy:
Excitation with evanescent field
- Advantages:
 1. better axial resolution
 2. better SNR, no fluorescence background
- Problem in optical design:
Extremely small illumination ring-shaped channel around the observation light cone



Epi-Fluorescence




TIRF



2 Optical System of the Microscopy II

Adjustment of Objective Lenses

- Adjustment of air gaps to optimize spherical aberration
- Reduced optimization setup



$$c_j = c_{j0} + \sum_{k=1,4} \Delta c_j \cdot \frac{\partial c_j}{\partial t_k}, \quad j=2,4,6,8$$


- Compensates residual aberrations due to tolerances (radii, thicknesses, refractive indices)

	d_2	d_4	d_6	d_8	C_{20}	C_{40}	C_{60}	C_{80}	W_{rms}
nominal	0.77300	0.17000	3.2200	2.0500	0.00527	-0.0718	0.00232	0.01290	0.0324
d_2 varied	0.77320	0.17000	3.2200	2.0500	0.04144	-0.07586	0.00277	0.12854	
d_4 varied	0.77300	0.17050	3.2200	2.0500	0.03003	-0.07461	0.00264	0.01286	
d_6 varied	0.77300	0.17000	3.2250	2.0500	0.00728	-0.07367	0.00275	0.01284	
d_8 varied	0.77300	0.17000	3.2200	2.0550	0.005551	-0.0717	0.00235	0.01290	
optimized	0.77297	0.16942	3.12670	3.2110	0.000414	0.00046	0.00030	0.01390	0.00468

2 Optical System of the Microscopy II

Adjustment of Objective Lenses


- Significant improvement for one wavelength on axis
- Possible decreased performance in the field



2 Optical System of the Microscopy II

Adjustment and Compensation

- Example microscopic lens
- Adjusting:
 1. Axial shifting lens : focus
 2. Clocking: astigmatism
 3. Lateral shifting lens: coma
- Ideal : Strehl $D_S = 99.62 \%$
With tolerances : $D_S = 0.1 \%$
After adjusting : $D_S = 99.3 \%$



2 Optical System of the Microscopy II Adjustment and Compensation

- Successive steps of improvements

FSF
(intensity normalized)

FSF
(energy normalized)

With Beams

Sp1 (Z_4 Z_9)

Sp2 (Z_7 Z_8)

